



TECNO TREND

Aseptic Powder Machines

| MICRO POWDER DOSING MACHINES (WITH ROTATING DISK) | | | |
|---|-------------------------------------|-------------------------------------|-------------------------------------|
| Technical data | FV100T | FV200T | FV300T |
| Output – single dose (vials/hr) | 6.000 | 12.000 | 18.000 |
| Vial dimensions diam. | 16 mm - 53 mm | 16 mm - 32 mm | 14 mm - 46 mm |
| Vial dimensions h | 40 mm - 110 mm | 2,5 times the diameter | 35 mm - 140 mm |
| Minimum dose (Min. dosages can change according to the powder specific weight) | 0,2 cc | 0,2 cc | 0,2 cc |
| Maximum dose. (the maximum volume is theoretically unlimited since it depends on dosages number, vial internal diameter and powder specific weight) | 7,5 cc (1,5 g. with single dosage) | 7,5 cc (1,5 g. with single dosage) | 7,5 cc (1,5 g. with single dosage) |
| Power consumption | 7 kw | 7 kw | 7 kw |
| Compressed air consumption | 50 Nl/min! | 50 Nl/min! | 50 Nl/min! |
| Vacuum pump | 100m ³ /h | 100m ³ /h! | 100m ³ /h! |
| Compressed air pressure in bar | 6 | 6 | 6 |
| PLC standard | SIEMENS | SIEMENS | SIEMENS |
| Net weight (kg) | 1100 | 1200 | 1300 |
| Gross weight (kg) | 1400 | 1500 | 1600 |